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Model CA-200+VA-200

Mitsubishi Chemical Corporation
Tokyo, Japan


LAB Online Exhibition

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
Model CA-200


Mitsubishi Chemical Corporation
Tokyo, Japan

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


No	名稱	功能
1	電源開關	By pressing [] side, power is on. By pressing [○] side, power is off.
2	5.7吋 LCD 螢幕顯示	It is 5.7 inches color LCD and displays setting and operation information and error messages.
3	操作面板 按鍵	For titration start, stop, and sample weight input. For keypad operation, display, and each input setting, refer to 5. Moisture Measurement Mode.
4	列印器	Titration parameters, measurement result, and graphs are printed.

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

No	名稱	功能
1	轉速調整鈕	For adjusting stirrer bar rotation speed Usually set it to "2 to 3". By setting it to "0", a stirrer bar stops
2	攪拌磁石	For stirring evenly anode solution (generator solution), samples, and iodine
3	樣品注入孔	For injecting liquid sample from here with a syringe
4	滴定杯	For Karl Fischer reaction
5	陽極電極 (檢測電極)	For potential detection
6	盲栓	Open a sealing plug and inject powdery and solid samples. Moisture can be prevented from pouring into a generator solution cell.
7	乾燥管盲栓	With a gap to discharge gas at much sample injection.
8	乾燥管(大)	For preventing moisture from pouring into a generator solution cell
9	陰極電極盲栓	For preventing moisture from pouring into a cathode solution cell
10	陰極電極槽 (雙白網電極)	For iodine generation from anode platinum mesh
11	攪拌器編號	Attach a stirrer number sticker.

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按 鍵	名 稱	功 能
	參數鍵/ 數字切換鍵	<ol style="list-style-type: none"> 1. For parameter setting Parameters can be set regularly and changed during titration. 2. The name input window is indicated at name input.
	樣品輸入鍵	<p>For sample parameter setting It can be set regularly. If required mass is not set, input is waited after titration.</p>
	功能鍵	<ol style="list-style-type: none"> 1. For each function selection 2. Unavailable during measurement.
	空白平衡 滴定鍵	<ol style="list-style-type: none"> 1. For setting titration to ON or OFF 2. For suspending titration and returning to the initial titration window
	開始/滴定鍵	<ol style="list-style-type: none"> 1. For starting moisture measurement and ending it forcibly 2. For canceling operations during dehydration waiting, delay, purge, pre-baking, and cooing to proceed to the next step 3. For stopping KF-200 piston up and down (at buret manual operation)
	列印鍵	<ol style="list-style-type: none"> 1. For printing statistical calculation result, file contents, and data 2. Unavailable during titration

按 鍵	名 稱	功 能
	上下移動鍵	<ol style="list-style-type: none"> 1. For transferring to items and selecting characters 2. For KF-200 piston up and down (at buret manual operation)
	左右移動鍵	<ol style="list-style-type: none"> 1. For changing parameters and selecting characters 2. For KF-200 valve selection (at buret manual operation)
	數字鍵	<ol style="list-style-type: none"> 1. For menu selection, titration and sample parameter input 2. For sample size decimal points
	符號鍵	For name input.
	跳出鍵	<ol style="list-style-type: none"> 1. For suspending the setting. 2. For returning to the previous window
	清除鍵	For deleting the right indication from the cursor position
	執行/確認鍵	For determining contents and menus
	頻道切換鍵	Available when a 2CH board is installed In channel selection, it beeps once for CH1 and it beeps twice for CH2.
	調查協助鍵	The help window is displayed to search troubleshooting information.
	列印紙輸出鍵	For paper feed Unavailable during printing

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Format 4

Format 5

Concentration calculation formulas

No.	Formula	Unit	Remark
0	M	H ₂ O μ g	Only moisture content
1	M/ (W-w)	Unit selection	Use it for measuring weighted liquid and solid samples.
2	M*X/ (W-w)		Use it when solving liquid or solid samples in solvents and measuring some samples. X: Dilution multiple
3	M/ (L*S)		Use it for weighing liquid samples by volume and measuring them L: Volume(ml), S: Density
4	M*Y*(1+t/ 273)/ G		Use it for weighing gaseous samples by volume and measuring them. t : Temperature (°C) G : Volume (l) Y : Input coefficient 1.244.
5	M-B		H ₂ O μ g
6	(M-B)/ (W-w)	Unit selection	Formula to subtract blank values from formula 1
7	(M-B)*X/ (W-w)		Formula to subtract blank values from formula 2
8	(M-B)/ (L*S)		Formula to subtract blank values from formula 3
9	(M-B)*Y*(1+t/ 273)/ G		Formula to subtract blank values from formula 4
10	Option Formula 1		For preparing formulas except 0 to 9 fixed ones
11	Option Formula 2		

Result units list

Indications	Contents
0: %, ppm Auto Select	Less than 0.1%: ppm is printed. 0.1% or more: % is printed.
1: % only	Only % is printed.
2: ppm only	Only ppm is printed.
3: Not Print	No unit print
4: Option Unit 1	Except the above units, an optional unit is prepared.
5: Option Unit 2	

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