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
OxiTop® IS 6 OxiTop® IS 12

BOD = Biochemical Oxygen Demand

Please read this information carefully before putting the device into service!

This device has been built and tested according to the IEC 1010 safety standards for electronic measuring instruments and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument can only be ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present operating instruction.

- Before connecting the plug-in power supply unit to the electricity supply network, it must be ensured that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).
- Attention. Magnetism! Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac, cardiac pacemakers...).
- The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of these operating instructions.
- When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.
- Balancing, maintenance and repair work must only be carried out by a suitably qualified technician authorized by us.
- If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.
- The safety of the user may be endangered, e.g., if the instrument
 - shows visible damage,
 - no longer operates as specified,
 - has been stored over a longer period under unsuitable conditions,
 - has been subjected to difficult conditions during transport.
- If in doubt, the instrument should as a rule be sent back to the manufacturer for repair and maintenance.

Contents

Instruction manual Manometric BOC Measuring Devices

OxiTop® IS 6 / IS 12

Recommended range of application	26
Delivery equipment.....	26

Instruction manual Operation of the Single Measuring

System OxiTop® 27

Measuring principle.....	27
BOD ₅ determination	28
Required instruments and accessories	28
Selection of sample volume	28
Sample preparation	28
Measurement.....	29
Disturbances.....	31
System messages	31
Cleaning	32
Battery replacement	32
Information.....	33
Technical data OxiTop® single measuring system.....	33

Instruction manual

Inductive Stirring System IS 6 / IS 12 34

Mains power supply	35
Operation notes	36
Information.....	36
Recommended accessories	37
Spare parts	37
Technical data	38
Test certificate EMC	39

Recommended range of application


OxiTop® measuring systems with inductive stirring system, for respirometric determination of the biochemical oxygen demand (BOD) according to the self-checking method.

Delivery equipment

Model	OxiTop® IS 6	OxiTop® IS 12
Order Nr.	208 210	208 211
Inductive stirring system (Stirring device)	6 OxiTop® instruments	12 OxiTop® instruments
Equipment	3 OxiTop® yellow 3 OxiTop® green	6 OxiTop® yellow 3 OxiTop® green
<u>Accessories:</u> Accessory set 1 (brown sample bottle, stirring rods, rubber quivers)	6 x	12 x
Accessory set 2 (stirring rod remover, overflow measuring beaker 164 ml and 432 ml, sodium hydroxide tablets, note block)	1 x	1 x

Operation of the Single Measuring System

OxiTop®

Measuring principle

BOD measurement with the OxiTop® measuring system is based on pressure measurement (difference measurement). The measuring is made by pressure measurement via **piezoresistive electronic pressure sensors**.

With the following functions, the OxiTop® measuring system minimizes the measuring work and is especially suited to the courses of the respirometric BOD measurement:

- **AUTO TEMP function:**

Control of the **temperature adjustment** and **automatic start** of the measurement (1 hour after starting the OxiTop® measuring system at the earliest, after 3 hours at the latest). It is not necessary to temper the samples to exactly 20°C before starting OxiTop®. Samples between 15°C and 20°C can be started immediately as BOD measurement will not begin until a stable 20°C is reached. This is controlled by the "AutoTemp" feature.

- **Data logging**

Automatic measurement with daily storage of measured value for 5 days: facilitates measurement without supervising, e. g. over the weekend.

- **Current value**

Indication of measured value (0 - 40 digits) corresponds to the digits of the previous WTW BOD instruments. Conversion of mbar to digits is stored in the instrument, i. e. evaluation of the measured values remains unchanged.

- **Measuring range reserve**

from 40 to 50 digits - for overranging without reset by opening the bottle.

BOD₅ determination

Normally municipal wastewater does not contain toxic or impeding substances. There are enough nutrient salts and suitable microorganisms. Under these conditions the **BOD₅ determination** with the OxiTop® measuring system is possible **in the undiluted sample**

Required instruments and accessories

- OxiTop® measuring system
- Inductive stirring system
- Incubator thermostatic box (temperature 20°C ± 1K)
- Sample bottles brown (nominal volume 510 ml)
- Stirring rods
- Stirring rod remover
- suitable overflow measuring beakers
- Rubber quivers
- Sodium hydroxide tablets

Selection of sample volume

Estimate the BOD₅ value to be expected for the wastewater sample:

Expected BOD₅ value ≈ 80 % of the COD value

Look for corresponding measuring range in the following chart and gather correct values for sample volume and factor.

Sample volume (ml)	Measuring range (mg/l)	Factor
432	0 - 40	1
365	0 - 80	2
250	0 - 200	5
164	0 - 400	10
97	0 - 800	20
43.5	0 - 2000	50
22.7	0 - 4000	100

Sample preparation


See following WTW application reports:

- WTW Application report 895230:
"BOD measurement in household wastewater"
- WTW Application report 895231:
"BOD measurement in organically heavily contaminated wastewater"
- WTW Application report 895232:
"BOD measurement with presence of impeding or toxic substances"

Measurement


Information: To measure the sample volume, usually overflow measuring beakers or measuring cylinders are used. Chose the volume according to the measured value expected. Too large measuring ranges will lead to inaccurate results. To estimate the measured value you can reckon with approx. 80% of the CBS value.

- For **sample preparation and filling of the measuring bottles** see the following instruction: DIN 38409 part 52: "Measurement of the oxygen consumption" as well as WTW application reports, chapter "Sample preparation and measurement".
- **Rinse** measuring bottle with **sample**. Empty thoroughly.
- **Exactly measure** the required oxygen-saturated (thoroughly homogenized) **quantity of the sample** according to information.
- Put the **magnetic stirring rod** into the bottle.
- Insert a **rubber quiver** in the neck of the bottle.
- Put **2 sodium hydroxide tablets** into the rubber quiver with a tweezers. (Caution: The tablets must never come into the sample!)
- Screw OxiTop® directly on sample bottle (tightly close).

Start measurement:

Press S and M simultaneously. (2 seconds) until the display shows 00.

Display: Stored values are deleted.

- Keep the measuring bottle with the OxiTop put on **for 5 days at 20°C** (e.g. in an incubator). After the meas. temperature has been reached (after 1 hour at the earliest, after 3 hours at the latest; **AutoTemp function**), the OxiTop automatically starts the measurement of the oxygen consumption .
- During the 5 days the sample is continuously stirred. The OxiTop **automatically stores one value every 24 hours** for 5 days. To have the current value shown press the M key.

Display current measured value:

Press M until measured value is displayed (1 second).

e.g.

- **Readout of the stored values** after the 5 days have passed.


Recall stored values:

Press S until measured value is displayed (1 second).

Scroll to next day by repressing the S key while the measured value is displayed (5 sec). Fast scrolling by repeatedly pressing the S key.

Measuring day either stored value, e. g.

or system message "memory empty" (F=no measured value):

- Convert the displayed measured value (digits) into the BOD value with the following table (Digits x Factor = BOD5 in mg/l):

Sample volume (ml)	Measuring range (mg/l)	Factor
432	0 - 40	1
365	0 - 80	2
250	0 - 200	5
164	0 - 400	10
97	0 - 800	20
43.5	0 - 2000	50
22.7	0 - 4000	100

Disturbances

- **Measured value remains below measuring range**
The display shows zero or too low a value.
The measuring equipment is not water-tight.
Check rubber quiver, screw top and bottle.
Insufficient sample pretreatment or preservation.
The temperature of the sample had not sufficiently been adjusted (< 15°C).
- **Measuring range exceeded.**
The measuring range chosen is too small. With very high values (> 2000 mg/l) we recommend to predilute the sample.
Nitrification inhibitor (allylthiourea) is missing or lacking.

Errors due to procedure have not been mentioned.

System messages

Memory empty (IF = measured value of day 1 is missing)

Change batteries (approx. every 3 years)

Value remains below measuring range (< 0 digits)

Value exceeds measuring range (> 50 digits)

Cleaning


Cleaning of sample bottles

- **Do not use disinfectants!**
(Disinfectants will kill the required microorganisms!)
- Remove **gross contaminations** mechanically, e.g. with a **brush**.
- Rinse the bottles with **clear water** or with **water of the next sample**.
(After using detergents rinse thoroughly! Detergent residues may disturb the BOD₅ determination!)

Cleaning of the OxiTop® Single Measuring System

- **Do not use alcohol or acetone!**
- Clean with a soft cloth and aqueous soap solution

Battery replacement (approx. every 3 years)

- 1 - Push in spring catch.
- Open upper lid.
- 2 - Remove batteries.
- Put in new batteries, pay attention to polarity!
(e.g. WTW-Modell **Batt/Oxitop**®, Order No. 209 012).
- Place upper lid with tap on locking pin (lower lid). Attention! Do not squeeze cable connection.
- Close upper lid (make the hook click).

Disposal note:

Please dispose of the lithium batteries at a battery collection point.

Information

Never use sealing grease or other lubricants for your OxiTop® instruments. Some of these products contain solvents which cause severe damages on the plastic housings. The tightness of the BOD bottles is completely sufficient without grease. Rough contaminations and particles on the contact surfaces of the rubber quiver and the OxiTop® must always be wiped off. WTW grants no guarantee for damages caused by sealing grease.

Labelling:

To label the bottles clearly, WTW offers a marking set (6 marking rings with numbers, which can be put on the neck of the bottle):

Model **MARK-6** (6 pcs)
Order No. 209 013.

Check of the measuring system:

On request WTW offer special test instruments for test control.

Further literature:


WTW BOD handbook
DIN 38409T51
DIN 38409T52

Technical data OxiTop® single measuring system

Measuring principle	Respirometric (manometric)
Display	2-digit 7-segment LED, 10 mm
Measured parameter	Digits
Measuring range	0 ... 40 digits (+10 digits for overranging)
Admissible measuring temperature	20°C ± 1 K
Admissible sample temperature (when being filled)	15-20°C
Power supply	Lithium batteries (260 mAh), 2 x type CR 2430 (3V)
Power consumption	Max. 25 mA (during measurement)
Protection class	3 IEC 1010
Protection system	IP 54 IEC 529
EMC Emissions	EN 50081-1, FCC Class A
Immunity	EN 50082-2, NAMUR recommendation
Climate class	2, VDI/VDE 3540
Ambient temperature	Storage: -25°C ... +65°C Operation: +20°C ± 1K
Test mark	CE
Dimensions	H: 69 mm, diameter 70 mm
Weight	Approx. 85 g

Inductive Stirring System

IS 6
IS 12

Mains power supply


Use original plug-in power supplies only (see chapter TECHNICAL DATA).

Before the instrument is plugged in and switched on, check that the working voltage of the instrument and line voltage are corresponding!

Connect power-supply unit with electric circuit and the instrument.

Operation notes

- Place instrument into incubator or thermostatic box.
- Connect line adaptor to mains socket.
- Connect line adaptor to instrument. After a short waiting time the Inductive Rotary field is automatically built up and the stirring rods start running.
- Place prepared sample bottles onto IS 6 / IS 12 inductive stirring system.

Information

- Interval-controlled Inductive Stirring System with program-controlled stirring procedure.
This means an optimum exchange between air space and sample takes place. The automatic, cycle-controlled stirring rod centering with synchronization of the rotary fields guarantees a safe stirring operation. Magnetic stirring rods fallen out of step are resynchronized. The sample bottles can be placed on IS6 / IS12 stirring systems even when they are in operation. The instruments operate without mechanically moved parts and are resistant to wear and maintenance-free.
- If the stirring operation does not start, disconnect and then reconnect the line adaptor to the mains socket. This will cause the stirring rods to move.
- **To ensure an optimum and secure stirring function use only original WTW stirring rods RST 600.**
- Cleaning with soft cloth and aqueous soap solution (**no alcohol or acetone!**).

Recommended accessories

	Model	Order-No.
Storing rack for 6 OxiTop® measuring systems	St-OxiTop®	209 010
Marking ring set (6 pcs.) with label set (1-12) for bottles	Mark-6	209 013
Spare batteries (1 set)	Batt/OxiTop®	209 012
Upgrade kit for 2 meas. places, with 1 of each OxiTop® yellow or green, 1 bottle, 1 small rubber bung, 1 stirring rod	OxiTop® SET-2	208 819
Overflow measuring beaker	Content 22.7 ml 43.5 ml 97.0 ml 164 ml 250 ml 365 ml 432 ml	MK 22/60 209 030 MK 43/600 209 040 MK 97/600 209 050 MK 164/600 209 060 MK 250/600 209 070 MK 365/600 209 080 MK 432/600 209 090

Spare parts

	Model	Order-No.
Spare OxiTop® measuring system yellow	OxiTop®/GB	208 810
Spare OxiTop® measuring system green	OxiTop®/GN	208 800
Sample bottle brown	PF 600	209 100
Stirring rods	RST 600	209 120
Rubber quiver	GK 600	209 170
Stirring rod remover	REF 600	209 130
Sodium hydroxide tablets	NHP 600	209 140
Nitrification restraint (ATH solution 250 ml)	ATH 600	209 330
Note block	Dia/600	209 020

Technical data

Power supply	Plug-in power supply unit (Euro, US, UK, Australian plug): Friwo FW7556 / 24 Input: 100 ... 240 V ~ / 50 ... 60 Hz / 400 mA Output: 24V DC / 750 mA
Power consumption	24 VA max.
Speed of the stirring rods	programm-controlled 180 ... 450 min ⁻¹
Instrument protection class	<ul style="list-style-type: none">• protection class 3, EN 61010-1• protection class IP30
EMC	<ul style="list-style-type: none">• Emissions EN50081-1/FCC class A• Immunity EN50082-1, EN50082-2, Namur
Ambient temperature	Storage: - 25°C ... + 65°C Operation: + 5°C ... + 40°C
Relative humidity	Yearly mean: < 75 % 30 days/year: < 95 % Remaining days: < 85 %
Test mark	CE
Dimensions (mm)	IS 6, IS 6 / 120 (B)265 x (T)181 x (H)67 IS 12, IS 12 / 120 (B)350 x (T)266 x (H)67
Weight (inclusive line adaptor)	IS 6, IS 6 / 120 1.6 kg IS 12, IS 12 / 120 3.2 kg

Test certificate EMC: IS 6

SENTON

EMV-Prüfzentrum ■ EMI/EMC-Testcenter ■

CERTIFICATE
to the Electromagnetic Compatibility
 to Test Report No. 52501-51217-1/-2

Sample: IS 6
Description: Stirrer
Applicant: WTW GmbH
Regulations: EN 50081-1:1992
 FCC part 15 subpart B
 EN 50082-2:1995
 Namur Recommendation 5/1993

Testresult:
 The samples are in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

RFI Emissions		
EN 55011:1991	Limit class B, group 1	
FCC part 15 subpart B	Limit class A	
Immunity		
Electrostatic Discharge IEC 801-2:1991	Contact discharge	4 kV
	Air discharge	8 kV
Electromagnetic Fields ENV 50140:1993	80 MHz - 1000 MHz	10 V/m
	Broadcasting frequency ranges	3 V/m
ENV 50204:1995	900 ± 5 MHz (pulse modulated)	10 V/m
Electrical Fast Transients IEC 801-4:1988	Power line (AC)	2 kV
Surge Immunity IEC 801-5	Power line - symmetrical coupling	1 kV
	Power line - unsymmetrical coupling	2 kV
Conducted RF immunity ENV 50141:1993	150 kHz - 80 MHz	10 Veff
	Broadcasting frequency ranges	3 Veff
Limitation of Starting Current Namur Rec. 5/1993	Starting current (peak)	≤ 15 x I _N
Mains Interruption Namur Rec. 5/1993	Interruption time	≥ 20 ms
Mains Supply Tolerance Namur Rec. 5/1993	Voltage	UN +20 %/-15 %
	Frequency	FN ± 6 %
50 Hz Magnetic fields EN 61000-4-8		30 A/m

Senton GmbH
 Johann Roidt

Straubing, 17. January 1996

Test certificate EMC: IS 12

SENTON

■ EMV-Prüfzentrum ■ EMI/EMC-Testcenter ■

CERTIFICATE
to the Electromagnetic Compatibility
 to Test Report No. 52501-51218-1/-2

Sample: IS 12
Description: Stirrer
Applicant: WTW GmbH
Regulations: EN 50081-1:1992
 FCC part 15 subpart B
 EN 50082-2:1995
 Namur Recommendation 5/1993

Testresult:

The samples are in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

RFI Emissions

EN 55011:1991	Limit class B, group 1	
FCC part 15 subpart B	Limit class A	
Immunity	(see test report no. 52501-51217-1 for IS 6)	
Electrostatic Discharge	Contact discharge	4 kV
IEC 801-2:1991	Air discharge	8 kV
Electromagnetic Fields	80 MHz - 1000 MHz	10 V/m
ENV 50140:1993	Broadcasting frequency ranges	3 V/m
ENV 50204:1995	900 ± 5 MHz (pulse modulated)	10 V/m
Electrical Fast Transients	Power line (AC)	2 kV
IEC 801-4:1988		
Surge Immunity	Power line - symmetrical coupling	1 kV
IEC 801-5	Power line - unsymmetrical coupling	2 kV
Conducted RF immunity	150 kHz - 80 MHz	10 Veff
ENV 50141:1993	Broadcasting frequency ranges	3 Veff
Limitation of Starting Current	Starting current (peak)	≤ 15 x I _N
Namur Rec. 5/1993		
Mains Interruption	Interruption time	≥ 20 ms
Namur Rec. 5/1993		
Mains Supply Tolerance	Voltage	UN +20 %/-15 %
Namur Rec. 5/1993	Frequency	FN ± 6 %
50 Hz Magnetic fields		30 A/m
EN 61000-4-8		

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What can Xylem do for you?

We're a global team unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to xylem.com.

Service and Returns:

Xylem Analytics Germany
Sales GmbH & Co. KG
WTW
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

Tel.: +49 881 183-325
Fax: +49 881 183-414
E-Mail wtw.rma@xylem.com
Internet: www.WTW.com

Xylem Analytics Germany GmbH
Dr.-Karl-Slevogt-Str. 1
82362 Weilheim
Germany

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SENTON GmbH · EMV-Prüfzentrum · Äußere Frühlingstraße 45 · 94315 Straubing · Tel. 00421/42081