LAB Online Exhibition

Melting Point Apparatus SMP3


Instructions for use

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English


Figure 1

Thank you for purchasing this Barloworld Scientific product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

Electrical Installation

THIS EQUIPMENT MUST BE EARTHED

Before connection please ensure that the mains power supply corresponds to that shown on the rating plate. Power consumption is 60W.

There is an IEC socket at the rear of the instrument for connection to the mains supply, (see figure 1). The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a UK 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains supply. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable you should obtain a moulded lead locally. If this is not possible, take the lead with the UK plug and replace the plug with a suitable alternative. See the enclosed instruction sheet for advice on how to carry out this procedure.

Should the mains lead require replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC 320 plug should be used.

N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

The mains lead should be connected to the instrument **BEFORE** connection to the mains supply.

**CAUTION: Fuse fitted in both live and neutral lines
IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

General Description

WARNING: The sample block may be hot.

The apparatus is designed to measure and record the temperatures of crystalline sample held within capillary tubes. Up to three tubes are accommodated in an illuminated chamber within an aluminium block. The tubes are viewed through a magnifying lens on the front of the unit. All controls are via a membrane keypad.

The temperature range is ambient to 360°C and the heating rate is variable between 0.5°C / minute and 10°C / minute. Temperatures are displayed on an LCD. Up to eight temperatures per capillary tube can be recorded and recalled on the LCD. The plateau temperature and heating rate are fully programmable. The SMP3 prompts can appear in English, French, Italian, Spanish and German.

When the melt cycle is started the SMP3 heats up to the pre-programmed plateau temperature at maximum rate. An alarm sounds and the SMP3 remains steady at the plateau temperature. With a touch of the start button, the SMP3 begins to heat at the pre-programmed rate.

Once the melt has occurred and the cycle stopped, cooling is automatic and can be selected to be either to the plateau temperature or to ambient.

Controls

Figure 1 shows the features for the SMP3.

Mains On/Off: This is a rocker type switch. Pressing the bottom switches the unit On. Pressing the top switches the unit OFF.

Sample Block: This aluminium block can hold three capillary sample tubes.

Viewer: Comprises of a magnifier lens encased in a plastic eye viewer.

Mains Indicator: This light illuminates green when the unit is on.

LCD: This display shows programming information and temperature readings.

Keypad: This membrane keypad controls the SMP3 via the following buttons:-

“MODE” key: Changes the display to select the plateau temperature and ramping rate.

“START” key: Starts the melting cycle.

“STOP” key: Stops the melting cycle.

▲“UP” key: Used in conjunction with the «MODE» key to select the plateau temperature and ramping rate. Also runs through the stored temperature readings in conjunction with the «STORE SAMPLE» keys.

▼“DOWN” key: Used in conjunction with the «MODE» key to select the plateau temperature and ramping rate.

“STORE SAMPLE 1 - 3” keys: Used to store and recall up to eight temperature readings per capillary tube.

Safety Advice Before Use

- ❖ The ventilation slots in the front of the unit should not be sealed or obstructed.
- ❖ Never move or carry the unit while in operation or when connected to the mains electricity supply.
- ❖ The unit should be carried using both hands with fingers under each side frame.

Operation

Make sure the mains lead is firmly connected in the IEC socket located at the rear of the unit and connect to the electricity supply. Switch the unit ON at the mains On/Off switch on the back of the instrument. When switched ON, a green light on the front panel and the sample chamber will illuminate.

The SMP3 prompts are multi language. To select a different language to English; hold down the «▲» and «▼» keys simultaneously then switch the unit ON at the mains. Select either ENGLISH, ESPANOL, FRANCAIS, DEUTSCH or ITALIANO using the «▲» key. When the required language is indicated press «START».


The display will first show “STUART SMP3” followed by “VERSION 3.0”. After a short pause, the display will alternate between;

The temperature displayed is the current ambient. The cooling pump can be heard. The plateau temperature and post plateau ramping rate may now be set.


Selecting the plateau temperature & ramping rate

Press the «MODE» key. The display will change to:

Use the «▲» and «▼» keys to select the desired plateau temperature (usually about 10°C below the approximate melting point).

If «MODE» is pressed again the screen will change to:

Use the «▲» and «▼» keys to select the required ramping rate. Press «MODE» once more to return to the “COOLING TO AMBIENT” screen.

Melting point measurement

1. Place a small amount of sample in a capillary tube and insert into the aluminium sample block. Up to three tubes can be accommodated at any one time.
2. Adjust the angle of viewing for maximum comfort.
3. Select the required plateau temperature and heating rate as described above.

- To begin the melt cycle press the «START» key. The pump will stop and the unit will heat at maximum rate to the pre-programmed plateau temperature. The screen will show "HEATING" along with the current temperature of the block. When the plateau temperature is reached an alarm will sound and the SMP3 will maintain that temperature. The display changes to "STEADY" followed by the current temperature.
- To continue heating at the pre-programmed ramping rate, press «START». The display will change to;



The temperature displayed is the current. The block is now being heated at the pre-programmed ramp rate.

N.B. The maximum temperature the block will reach in this phase is 10°C above the plateau temperature plus 20 times the programmed rate.

- While the SMP3 is in the ramping mode, temperatures can be recorded via the «STORE SAMPLE» keys (see below).

When the melt has occurred press «STOP». The display will change to "COOLING TO AMBIENT" and the cooling pump can be heard. Now the stored temperatures can be viewed (see below).

- If another test at the current plateau temperature and ramping rate is to be carried out, pressing «START» will tell the SMP3 to cool to the plateau temperature. The screen will alternate between;

Once the plateau temperature has been reached the display will change to "STEADY" and the alarm will sound.

Begin ramping by pressing «START».

- If another run with different parameters is required, select a new plateau temperature and ramping rate as described above. Press «START» to reach the new plateau temperature.

Recording temperatures


- Temperatures are recorded via the «STORE SAMPLE» keys. There are three memories (one for each tube) and each memory can hold eight temperature values.

N.B. If a ninth value is recorded an alarm will sound to indicate that the memory is full.

- Temperature values can be recorded in the "RAMPING" mode by pressing the appropriate «STORE SAMPLE» key.

Viewing stored temperatures

- The stored value(s) can be viewed during the "COOLING TO AMBIENT" mode. When one of the «STORE SAMPLE» keys is pressed the screen will show;

Here «STORE SAMPLE 1» has been pressed and the first temperature inputted was 202.1°C. Pressing the «▲» key will show the other seven inputs in turn. Pressing «STOP» will return to the "COOLING TO AMBIENT" mode.

N.B. Temperature values remain in the memory until another melt cycle is initiated by pressing «START».

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

The magnifier lens can be removed from its retainer and gently cleaned with a damp cloth. To remove the test block window, slip out the window retaining circlip. Remove the window using tweezers and clean with a damp cloth.

Any repairs or replacement of parts **MUST** be undertaken by **suitably qualified personnel**.

Only spare parts supplied or specified by Barloworld Scientific or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, contact the Technical Department of Barloworld Scientific or the point of sale.

Spares and Accessories

For a comprehensive list of parts required by service engineers conducting internal repairs, please contact the Sales Department of Barloworld Scientific quoting both the model and serial number.

Technical Specification

Number of samples	Three
Temperature range	Ambient to 360°C
Temperature resolution	0.1°C
Temperature accuracy	± 0.5°C at 20°C ± 1.5°C at 300°C
Display	Alphanumeric LCD
Sensor	PT100
Ramp rate	0.5° to 10°C minute
Memory storage	8 x temperature readings per sample tube tested
Dimensions (wxdxh)	180 x 260 x 300 mm
Net weight	3.2 Kg

Warranty

Barloworld Scientific Ltd. warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of **three (3)** years. In the event of a justified claim Barloworld Scientific will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by incorrect installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards

and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Barloworld Scientific Ltd.

Declaration of Conformity

Melting Point Apparatus

Model SMP3

This product complies with the requirements of the EU Directives listed below:

89/336/EEC Electromagnetic Compatibility Directive, amended by 92/31/EEC and 93/68/EEC.

73/23/EEC Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326: 1997 + Amendments A1: 1998, EN61000-4-2: 1995, EN61000-4-4: 1995, EN61000-4-5: 1995, EN 61000-4-6: 1996, EN61000-4-9: 1993, EN61000-4-11: 1994, EN61000-3-2: 1995, EN61000-3-3: 1995, EN 61010-1: 1993 + Amendments A1: 1992 and A2: 1995, EN61010-2-010: 1994.

Compliance Certificates and Full Reports.

Ref: RETS0287/11 and RETS0287/12

From "Epsilon Technical Services" an Independent Accredited Test House Showing Compliance to the above Standards are available on request.

CE Mark Affixed 01

Signed:  (Mr D E Hicks)

Date: 25/7/05

Authority: Technical / Development Manager

for

Barloworld Scientific Limited - Stone - Staffs - ST15 0SA - UK
Tel: +44 (0) 1785 812121 Fax +44 (0) 1785 813748

INSPECTION REPORT

MODEL SMP3


ELECTRICAL SAFETY

1. Earth continuity 
2. Insulation 
3. Flash test 

FUNCTIONAL

1. Indicators 
2. Temperature control 
3. Visual acceptance 

QUALITY CONTROL INSPECTOR

**Barloworld
Scientific**

Barloworld Scientific France SAS

ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bibby@bibby-sterilin.fr

Barloworld Scientific Italia Srl

Via Alcide de Gasperi 56
20077 Riozzo di Cerro al Lambro
Milano Italia
Tel: +39 (0)2 98230679
Fax: +39 (0)2 98230211
e-mail: marketing@barloworld-scientific.it
www.barloworld-scientific.it

Barloworld Scientific Ltd

Beacon Road Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
www.barloworld-scientific.com

Barloworld Scientific US Ltd

350 Commerce Drive
Rochester NY 14623
USA
Tel: (800) 828-6595
Fax: (585) 334-0241
www.dynalabcorp.com

Afora S.A.

Calle Aribau 240
08006 Barcelona Spain
Tel: +343 93-306 98 00
Fax: +343 93-306 98 23
e-mail: marketing@afora.com
www.afora.com