



敝公司提供電熱包系列產品，被廣泛使用於實驗室及產業界，而其表面最高溫度達攝氏450度，並接受客戶特殊規格品製作。

Heating Mantle with Stirrer 電熱包攪拌器

MODELS	CAPACITY	POWER
MNS-250	250 ml	200 w
MNS-500	500 ml	300 w
MNS-1000	1000 ml	400 w
MNS-2000	2000 ml	500 w
MNS-5000	5000 ml	700 w
MNS-12LC	12 L	1300 w
★附速度控制(0~1200rpm)及比例式調溫器		



Mantle Heater For Beaker With Controller 附調溫器燒杯型電熱包

MODELS	CAPACITY	POWER
BMN-250C	250 ml	200 w
BMN-500C	500 ml	300 w
BMN-1000C	1000 ml	400 w
BMN-2000C	2000 ml	500 w
BMN-3000C	3000 ml	600 w
BMN-5000C	5000 ml	700 w



Mantle Heater For Flask 軟式電熱包

MODELS	CAPACITY	POWER
SMN-50	50 ml	60 w
SMN-100	100 ml	80 w
SMN-250	250 ml	200 w
SMN-500	500 ml	300 w
SMN-1000	1000 ml	400 w
SMN-2000	2000 ml	500 w
SMN-3000	3000 ml	600 w
SMN-5000	5000 ml	700 w

+°C

Heating Mantles